

DECTRIS[®]

detecting the future

MYTHEN

PILATUS


EIGER

X-ray detectors for plasma diagnostics

Tilman Donath

DECTRIS Ltd.
5400 Baden
Switzerland
www.dectris.com

History of DECTRIS




2006: PILATUS 6M. First large-area hybrid-pixel X-ray detector

Facts:

- Spin off from PSI founded in 2006
- Located in Baden, Switzerland: 1500 m² production facilities, labs and office space (production, sales, and support)
- Products: 1D and 2D Hybrid-Pixel X-ray Detectors for scientific & industrial applications
- World leader in hybrid pixel X-ray detectors

- Staff by the end of 2013: ~50 FTE
- Turnover 2012: ~22 MCHF (17.5 MEUR)
- 151 systems delivered
- Strong focus on R&D! (ASIC design)
R&D expenses: 4.25 MCHF

Hybrid Pixel Technology



**PILATUS module
100'000 pixels**

Single-photon counting => Noise-free imaging

Large energy range: 1.8 – 30 keV

Radiation hard design

Business Segments

synchrotron


– **Highest performance**

laboratory & industry

– **High performance**

specific solutions

– **Special detectors**



PILATUS 2M at protein crystallography beamline



Mainly OEM partners. Example:
Rigaku - Compact HomeLab™
with PILATUS 100K




Customized PILATUS 900K
for EAST tokamak, IPP

Plasma spectroscopy systems with PILATUS

	System	Detector	Application
Alcator C-Mod (Boston, USA)	XICS Ar16+ system	3x Pilatus 100K	Routine analysis of plasma temperature and rotation
	XICS swappable crystal system	Pilatus 100K	Impurity density and high temperature studies <i>Both XICS systems operated in one Helium box</i>
	Energy Resolving Pinhole camera	Pilatus 100K	Demonstration system for impurity density and electron temperature imaging <i>Operated in air</i>
LHD (Japan)	XICS Ar16+ system	Pilatus 300K, water-cooled	Routine analysis of plasma temperature and rotation
	XICS Ar17+ system	Pilatus 100K camera	Routine analysis of plasma temperature and rotation in high temperature plasma <i>Both systems operated in vacuum</i>
EAST (China)	XICS Ar16+ system	2x Pilatus 100K cameras	Routine analysis of plasma temperature and rotation
	XICS multi crystals system	Pilatus 900K, specific in-vacuum detector	<i>delivered 07/2013</i>
KSTAR (Korea)	XICS Ar16+ system	Pilatus 100K	Routine analysis of plasma temperature and rotation
NSTX (Princeton, USA)	XICS Fe24+ system	2x Pilatus 100K	Impurity density, plasma temperature. (operated in 2010-2011, detectors were later shipped to EAST)

X-ray Imaging Crystal Spectrometer at Alcator C-Mod



- Bent-crystal spectrometer measuring Ar^{16+} emission line (3.1 keV)
- Operation in He-filled chamber
- Spatial resolution (1D poloidal profile) and energy resolution
- Plasma velocity (v) and temperature (T) are measured via the Doppler effect

HIREXSR system at Alcator C-Mod, 2007
Images courtesy of Matt Reinke

PILATUS3 900K-IPP for EAST tokamak



Detector:

- **9 modules**
(84 × 325 mm²)
- **in-vacuum**
operation
- **fast framing,**
up to 500 Hz
- **readout time**
0.95 ms/frame

delivered: 07/2013

commissioned:
09/2013

Description of XICS multi-crystal system design:
B. Lu *et al.*, Review of Scientific Instruments, **83**, 10E130 (2012)

DECTRIS[®]

detecting the future

***Thank you for
your attention!***

www.dectris.com

DECTRIS Ltd.
5400 Baden
Switzerland
www.dectris.com